

Body Assist Testicle Support - Suspensory

- For post-surgical conditions and before injury
- Ventilated elastic for cool, comfortable and extended use
- Helps to relieve the strain and fatigue of swollen testicles



- ✓ Contracted DVA product
- ✓ Contract approved

Code	Display Name
FBS303496	BodyAssist - Testicle Support - Suspensory - Xlarge
FBS303494	BodyAssist - Testicle Support - Suspensory - Small
FBS303495	BodyAssist - Testicle Support - Suspensory - Medium
FBS284420	BodyAssist - Testicle Support - Suspensory - Large