

## ROHO - Mini Max



- Designed for short-term needs of users with high comfort needs but at low-level risk of skin breakdown
- Inter-connected air cells provide shock absorption and pressure reduction
- Lightweight: At 0.9 kg the Mini Max cushion adds minimal weight to the entire mobility solution
- Great for travel and sports
- All models include: 2 way stretch cover, Hand pump & Repair kit, Operating instructions

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Manage tissue integrity

**Category:** Prescribed

Click [here](#) for more information

## Warranty Information

- 1 Year

permobil

**ROHO**

### NDIS

**Item Code:**

05\_181006190\_0105\_1\_2

**Pressure Care Cushion –**

**Composite:**05\_429\_0103\_1\_2

Click [here](#) for more information

### Compliance

EN-ISO 10993-10 (1999-11) / EN 1021-1 and EN 1021-2

| Code      | Display Name              | Total Length | Total Width | Total Height | Safe Working Load | Material  | Pressure Care Rating |
|-----------|---------------------------|--------------|-------------|--------------|-------------------|-----------|----------------------|
| PTC585800 | ROHO Mini Max - 15 x 15in | 15 in        | 15 in       | 1.25 in      | Unlimited         | Air Cells | Low Risk             |

| Code      | Display Name              | Total Length | Total Width | Total Height | Safe Working Load | Material  | Pressure Care Rating |
|-----------|---------------------------|--------------|-------------|--------------|-------------------|-----------|----------------------|
| PTC585900 | ROHO Mini Max - 16 x 16in | 16 in        | 16 in       | 1.25 in      | Unlimited         | Air Cells | Low Risk             |
| PTC586000 | ROHO Mini Max - 18 x 16in | 16 in        | 18 in       | 1.25 in      | Unlimited         | Air Cells | Low Risk             |