

M1 Corpus





- Less fuss, more adaptable: A quality solution for those users who require a more traditional rehab style seating system with limited to no power seat functions. Designed with adaptable seating and positioning needs in mind, the M1 offers all the flexibility without the compromise.
- M1 Core package available View Package Here
- Locking flip back arms with easy to release mechanism.
- Multi-angle adjustable swing-away legrests.
- FlexLink™ suspension and precision engineered CasterGlide™ technology to help ensure smooth, confident operation in a variety of driving scenarios.
- Designed with 1" diameter rigid tubular backrest frame allows mounting of virtually any seating and positioning solution.
- Transport friendly features including standard crash tested transport brackets, optional EZ fold backrest and a lower overall chair weight of 145 kg with batteries, expands available travel options

Warranty Information

• 1 year

\bigcirc	
Contract approve	ed

ndis

Item Code:

05_122306125_0105_1_2

Wheelchair Powered -

Composite:05_419_0105_1_2 Click here for more information

Code	Display Name	Seat Width	Total Width	Safe Working Load	Maximum Speed	Range
PWC652575	M1	14 - 22 in	24 in	136 kg	8 km/h	25 km