

Quickie AllCourt Ti Wheelchair



- The Quickie All Court® wheelchair is the new standard in performance and style for the recreational or elite basketball player
- Optimised for maneuverability and efficiency without adding complexity or extra parts
- The anti-tip brace locks out the rearward flex in the frame, giving you more power when you lean back to shoot over your competition

 Contract approved

Support @ Home

Function: Manual wheelchairs

Category: Prescribed

Click [here](#) for more information



Code	Display Name	Safe Working Load	Seat Depth	Seat Width
MWC427337	Quickie Allcourt Ti Wheelchair (As Per Script)	113 kg	255 - 510 mm	305 -510 mm