

Quickie AllCourt Ti Wheelchair



- The Quickie All Court® wheelchair is the new standard in performance and style for the recreational or elite basketball payer
- Optimised for maneuverability and efficiency without adding complexity or extra parts
- The anti-tip brace locks out the rearward flex in the frame, giving you more power when you lean back to shoot over your competition



Code	Display Name	Safe Working Load	Seat Depth	Seat Width
MWC427337	Quickie Allcourt Ti Wheelchair (As Per Script)	113 kg	255 - 510 mm	305 -510 mm