

ROHO - Single Compartment



permobil

ROHO

This cushion offers an optimal environment for individuals with skin/soft tissue integrity issues. The cellular design and the soft pliable material allow the individual to be immersed into the cushion for maximum pressure redistribution, while decreasing friction and shear. The simple adjustment allows the cushion to fit the shape of the individual, minimising tissue deformation and maximising blood flow.

- **Adjustable:** Adjusts to individual's shape and weight, and changes in individual's position over time
- **Lightweight:** Easy to handle and carry
- **Washable:** Cushion can be easily cleaned with soap and water. Cover is machine washable
- **Repairable:** The flame-resistant, black neoprene cushion can be patched or repaired to extend its use
- **Skin Protection:** ROHO DRY FLOATATION® Technology provides unparalleled performance in facilitating blood flow by fitting, matching and tracking the shape of a person
- **Convenience:** All postural adjustments are made while the user is seated, eliminating the need for transfer on and off the chair during the evaluation process
- Available in Low Profile (Mid and High Profile available as {Sensor Ready {link to NEWCODE1}})

Warranty Information

- 2 Years

- ✓ Contracted DVA product
- ✓ Contract approved

ndis

Item Code:

05_181006190_0105_1_2

Pressure Care Cushion –

Composite:05_429_0103_1_2

Click [here](#) for more information

Compliance

EN-ISO 10993-10 (1999-11) / EN 1021-1 and EN 1021-2

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Safe Working Load	Material	Pressure Care Rating
PTC544930	ROHO Single Compartment - LP - 11.5 x 11.5in (6 x 6cell)	1199 -		11.5 in	11.5 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC545000	ROHO Single Compartment - LP - 11.5 x 13.25in (6 x 7cell)	1199 -		13.25 in	11.5 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC545110	ROHO Single Compartment - LP - 13.25 x 13.25in (7 x 7cell)	1199 -		13.25 in	13.25 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC545130	ROHO Single Compartment - LP - 13.25 x 15in (7 x 8cell)	1199 -		15 in	13.25 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC545210	ROHO Single Compartment - LP - 13.25 x 16.75in (7 x 9cell)	1199 -		16.75 in	13.25 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC545600	ROHO Single Compartment - LP - 15 x 15in (8 x 8cell)	1199 -		15 in	15 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC545800	ROHO Single Compartment - LP - 15 x 16.75in (8 x 9cell)	1199 -		16.75 in	15 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC545300	ROHO Single Compartment - LP - 15 x 18.5in (8 x 10cell)	1199 -		18.5 in	15 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC546510	ROHO Single Compartment - LP - 16.75 x 15in (9 x 8cell)	1199 -		15 in	16.75 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC546700	ROHO Single Compartment - LP - 16.75 x 16.75in (9 x 9cell)	1219 -		16.75 in	16.75 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC546000	ROHO Single Compartment - LP - 16.75 x 18.5in (9 x 10cell)	1199 -		18.5 in	16.75 in	2.25 in	Unlimited	Air Cells	Very High Risk

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Safe Working Load	Material	Pressure Care Rating
PTC546200	ROHO Single Compartment - LP - 16.75 x 20in (9 x 11cell)	1199 -		20 in	16.75 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC543990	ROHO Single Compartment - LP - 18.5 x 15in (10 x 8cell)	1199 -		15 in	18.5 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC544100	ROHO Single Compartment - LP - 18.5 x 16.75in (10 x 9cell)	1219 -		16.75 in	18.5 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC543800	ROHO Single Compartment - LP - 18.5 x 18.5in (10 x 10cell)	1219 -		18.5 in	18.5 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC543910	ROHO Single Compartment - LP - 18.5 x 20.25in (10 x 11cell)	1199 -		20.25 in	18.5 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC544410	ROHO Single Compartment - LP - 20.25 x 16.75in (11 x 9cell)	1199 -		16.75 in	20.25 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC544350	ROHO Single Compartment - LP - 20.25 x 20.25in (11 x 11cell)	1199 -		20.25 in	20.25 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC544250	ROHO Single Compartment - LP - 20 x 18.5in (11 x 10cell)	1199 -		18.5 in	20 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC544550	ROHO Single Compartment - LP - 22 x 18.75in (12 x 10cell)	1369 -		18.75 in	22 in	2.25 in	Unlimited	Air Cells	Very High Risk
PTC544420	ROHO Single Compartment SR - HP - 22 x 16.75in (12 x 9cell)	1369 1		16.75 in	22 in	4.25 in	Unlimited	Air Cells	Very High Risk