

iCare Classic ActiveX Pillow



- The iCare Visco Classic pillow is made up of a high-density ActiveX™ core in a traditional classic pillow shape, a cooling gel insert, and a cotton cover pillow.
- Zip off washable cover
- Breathable vant fabric to increase airflow
- Suitable for pain management and pressure redistribution
- Suits front, back and side sleepers.
- Gel insert assists with temperature regulation

✔ Contract approved

Support @ Home

Function: Body stabilisation

Category: Under advice

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Total Height
BEP051055	iCare Classic ActiveX Pillow	600 mm	420 mm	110 mm