

## Aspire Overbed Table - Thermoform Recessed Top



- Recessed Thermoform top to prevent spills
- Lightweight mobile table positions overbed or chair
- Can be used for eating, drawing or other activities
- Flat top ideal for hospital or home use
- Fully height adjustable
- “C” base shape fitted with four castors enables close positioning to the bedside or chair
- Powder coated mild steel frame construction
- Cream standard laminate top with vinyl wrap around the edge, or “thermoform” cream top which is heat resistant

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Tables

**Category:** Prescribed

Click [here](#) for more information

### Warranty Information

- 1 Year

### NDIS

**Item Code:**

05\_180315322\_0103\_1\_2

**Minor Bed – Composite:**

05\_428\_0103\_1\_2

Click [here](#) for more information



### Compliance

ISO134185:2003 / EN ISO

13485:2012

Code	Display Name	Total Length	Total Width	Total Height	Material
BEA028000	Aspire Overbed Table - Thermoform Recessed Top - Cream	430 mm	800 mm	730 - 1095 mm	Thermoform Top

