

Teal Ottoman Box Bed

Robust construction, engineered for demanding mental health and challenging environments.

- Durable melamine faced chipboard construction (18mm thick)
- Fully sealed box bed solution
- Reinforced with security screws, ideal for extreme environment use
- Economical solution without compromising client and user safety



Warranty Information

- 1 Year

 Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
BEB050510	Teal Ottoman Box Bed - (Choose colour)	1985 mm	985 mm	370 mm	320 kg	70 kg