

## Teal Ottoman Box Bed



Robust construction, engineered for demanding mental health and challenging environments.



- Durable melamine faced chipboard construction (18mm thick)
- Fully sealed box bed solution
- Reinforced with security screws, ideal for extreme environment use
- Economical solution without compromising client and user safety

### Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
BEB050510	Teal Ottoman Box Bed - (Choose colour)	1985 mm	985 mm	370 mm	320 kg	70 kg