

Freedom Split Seat Chair - 3 in 1 - HBA370



- Patented split padded seat allows easy access to clean body parts from the front or rear of the chair without requiring the patient to stand
- Anti-slip suction rubber tips
- Adjustable legs with adjustment pins spaced 25 mm apart on three legs and 12.5mm apart on the fourth to ensure maximum stability. The oval shaped tubing prevents misalignment between the inner and outer tubes
- Powder coated aluminium frame
- Purchased separately:
 - Freedom Backrest - BTT163120 for added support
 - Freedom Splash Guard - GSP650720

✓ Contracted DVA product

✓ Contract approved

Compliance

AS/NZS 3973:2009

Warranty Information

- 1 Year



Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height
BTT163100	Freedom Split Seat Chair - 3 in 1 - HBA370	317	1	595 mm	645 mm	650 - 800 mm	160 kg	6.3 kg	450 mm	480 mm	440 - 590 mm
BTT163120	Freedom Backrest - HBA371	63	1	-	-	-	-	-	-	-	-