


Aspire MAXI LIFESTYLE Community Bed





- The Aspire Maxi Lifestyle Community Bed is a stylish, convenient and economical solution for users weighing up to 350 kg
- Small footprint when folded improves storage options in space constrained facilities
- Single button press electrical folding mechanism requires no lifting and improves workplace health & safety outcomes for delivery staff and maintenance teams
- **Transportable Design:** Rapid tool-free folding and unfolding for efficient setup/removal and transportation
- Small footprint: when folded improves storage options in space constrained facilities
- **Super-King size mattress platform:** caters to true bariatric clients reducing falls risk and moving and handling complexity.
- 2060 mm long sleeping surface with 2150 mm internal sleep area caters to taller patients
- Four Total-Lock Castors: Large 120 mm castors improve ease of movement compared to competing products
- Auto-regression in the backrest reduces skin shear forces and removes a significant cause of patient slippage. This reduces the need for care staff to reposition oversized patients
- Removable Folding Side Rails: purchased separately are durable enough to handle the rigors of equipment rental pools
- Accessories purchased separately <u>click here</u>
- Leading Quality Frame & Electronics: Linak electronics offer market leading reliability and performance
- Compatible with: <u>Aspire Bari10</u> alternating turning mattress (dynamic) and Aspire MAXI Lifestyle foam mattress (static) offering a high care rental package at a market leading price



Accessories

Click on the cart icon to display the range of optional accessories



- **Stylish Timber End Shrouds:** Are included as standard and offer a simple solution for cleaning and infection control
- **Highly Adjustable:** Four electric profiling functions to suit all comfort and positioning needs
- Full Safe Working Load: 385 kg, Maximum Patient Weight 350 kg

Warranty Information

Frame: 3 YearsElectronics: 2 YearsHandset: 1 year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Number of Functions
BEB046020	Aspire MAXI LIFESTYLE Community Bed	2290 mm	1360 mm	330 - 730 mm	385 kg	4