

Aspire MAXI LIFESTYLE Community Bed



- The Aspire Maxi Lifestyle Community Bed is a stylish, convenient and economical solution for users weighing up to 350 kg
- Small footprint when folded improves storage options in space constrained facilities
- Single button press electrical folding mechanism requires no lifting and improves workplace health & safety outcomes for delivery staff and maintenance teams
- **Transportable Design:** Rapid tool-free folding and unfolding for efficient setup/removal and transportation
- **Small footprint:** when folded improves storage options in space constrained facilities
- **Super-King size mattress platform:** caters to true bariatric clients reducing falls risk and moving and handling complexity.
- 2060 mm long sleeping surface with 2150 mm internal sleep area caters to taller patients
- **Four Total-Lock Castors:** Large 120 mm castors improve ease of movement compared to competing products
- **Auto-regression** in the backrest reduces skin shear forces and removes a significant cause of patient slippage. This reduces the need for care staff to reposition oversized patients
- **Removable Folding Side Rails:** purchased separately are durable enough to handle the rigors of equipment rental pools
- Accessories purchased separately [click here](#)
- **Leading Quality Frame & Electronics:** Linak electronics offer market leading reliability and performance
- **Compatible with:** [Aspire Bari10](#) alternating turning mattress (dynamic) and Aspire MAXI Lifestyle foam mattress (static) offering a high care rental package at a market leading price

 Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Accessories

Click on the cart icon to display the range of optional accessories

- **Stylish Timber End Shrouds:** Are included as standard and offer a simple solution for cleaning and infection control
- **Highly Adjustable:** Four electric profiling functions to suit all comfort and positioning needs
- Full Safe Working Load: 385 kg, Maximum Patient Weight 350 kg

Warranty Information

- Frame: 3 Years
- Electronics: 2 Years
- Handset: 1 year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Number of Functions
BEB046020	Aspire MAXI LIFESTYLE Community Bed	2290 mm	1360 mm	330 - 730 mm	385 kg	4