

iCare Companion Bed



- Suitable for use with the IC222 and IC333, the Companion Bed is perfect for partners, parents or caregivers
- The companion bed is a non-mechanical upholstered slat base, which sits beside an adjustable homecare bed

✔ Contracted DVA product

✔ Contract approved

Warranty Information

- 5 Years

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_181210200_0103_1_2

Power Bed – Composite:

05_423_0103_1_2

Click [here](#) for more information

Compliance

IEC 60601-2-52:2009

Code	Display Name	Total Length	Total Width	Safe Working Load
BEB045959	Companion Bed - Base Only - King Single	2150 mm	1120 mm	100 kg
BEB045961	Companion Bed - Base Only - Long Single	2150 mm	960 mm	100 kg

