

High Back Chair (Height adjustable)



- A high back orthapaedic chair that offers adjustability, comfort and support
- Contoured Backrest High back design provides increased postural support
- Softly padded armrests, seat and back upholstery are comfortable and supportive
- Easy clean vinyl surfaces
- Height Adjustable Easy push button adjustment for comfortable positioning and changing needs
- Sturdy Design Generous weight capacity, yet lightweight for easy manoeuvrability
- Non-Slip Tips Provide stability whilst sitting or transferring

*** Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

Code

Display Name

H_CHP208050 High Back Chair (Height adjustable)