

## High Back Chair (Height adjustable)



- A high back orthopaedic chair that offers adjustability, comfort and support
- Contoured Backrest – High back design provides increased postural support
- Softly padded armrests, seat and back upholstery are comfortable and supportive
- Easy clean vinyl surfaces
- Height Adjustable – Easy push button adjustment for comfortable positioning and changing needs
- Sturdy Design - Generous weight capacity, yet lightweight for easy manoeuvrability
- Non-Slip Tips – Provide stability whilst sitting or transferring



\*\*\* Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

**Code**

**Display Name**

H\_CHP208050 High Back Chair (Height adjustable)