

## Nordic Mobicline Highback Mobile Chair w/ Recline



- Manual recline mechanism
- Fully welded tubular frame
- Push handle and pull out footrest
- A large activity tray included
- 3 segment back cushion for comfort and individual fit
- Removable cushions for ease of cleaning
- Lockable 100 mm castors

 Contract approved

### Support @ Home

**Function:** Sitting furniture

**Category:** Under advice

Click [here](#) for more information

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Backrest Height
CHP194700	Nordic Mobicline Highback Mobile Chair w/ Recline - Vinyl Latte	690 mm	1120 mm	140 kg	490 mm	530 mm	680 mm
CHP194715	Nordic Mobicline Highback Mobile Chair w/ Recline - Ink Vinyl	690 mm	1120 mm	140 kg	490 mm	530 mm	680 mm