

Jobst Elvarex Armsleeve w/ Shoulder Strap



Jobst Elvarex Armwear Overview:

- Jobst Elvarex garments feature a special low-stretch, double-knit fabric, which exerts a light massaging effect on the initial lymphatics, and minimizes re-oedematisation
- The high cotton content, flat knit construction and short-stretch thread means that Elvarex products provide an advanced anatomical fit with maximum comfort
- Compression: Class 1: 20 - 30mmHg

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information



Code	Display Name	Calf/Arm Circumference	Elbow Circumference	Wrist Circumference
PTS642155	JOBST Elvarex - Armsleeve with shoulder strap (Right) - 20-30mmHg - Sml	230 - 260 mm	240 - 260 mm	160 - 180 mm
PTS642156	JOBST Elvarex - Armsleeve with shoulder strap (Left) - 20-30mmHg - Sml	230 - 260 mm	240 - 260 mm	160 - 180 mm
PTS642157	JOBST Elvarex - Armsleeve with shoulder strap (Right) - 20-30mmHg - Med	260 - 280 mm	260 - 290 mm	180 - 200 mm
PTS642158	JOBST Elvarex - Armsleeve with shoulder strap (Left) - 20-30mmHg - Med	260 - 280 mm	260 - 290 mm	180 - 200 mm
PTS642159	JOBST Elvarex - Armsleeve with shoulder strap (Right) - 20-30mmHg - Lge	280 - 300 mm	290 - 320 mm	200 - 230 mm
PTS642160	JOBST Elvarex - Armsleeve with shoulder strap (Left) - 20-30mmHg - Lge	280 - 300 mm	290 - 320 mm	200 - 230 mm