

Belmont Bedside Cabinet



- Sleek and stylish bedside cabinet
- Classic design will complement any decor
- Features an open shelf to keeps books, glasses and other personal items close at hand
- Lockable drawer keeps valuables and medications safe
- Small slimline handle features a curved shape that's easy to grab onto
- Narrow design suits smaller spaces
- Dimensions: 700 mm (H) x 405 mm (L) x 450 mm (W)
- Fabric colour/material to be advised at time of order



Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height
BEA828420	Belmont Bedside Cabinet	495mm	450mm	700mm
BEA826230	Belmont Bedside Cabinet - Timber	450 mm	470 mm	715 mm
BEA826360	Chelsea Bedside Cabinet - Laminate	450 mm	470 mm	715 mm