

Teal Roku Bedside Cabinet



Roku bedroom is a secure, extremely durable one-piece rotationally moulded collection for challenging environments.



- Smooth, high-impact polyethylene, continuously moulded for maximum durability
- Stain-resistant, seamless, and easy to clean for infection management
- Bedside cabinet sides fit flush to walls
- Smooth, round profile edges ensure enhanced safety and security
- 100% recyclable material

Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Unit Weight
BEA828500	Teal Roku Bedside Cabinet - (Choose colour)	460 mm	510 mm	595 mm	12.5 kg