

## Oppo 1080 Tennis Elbow Brace



- Provides full support of elbow region in conjunction with reducing of the insertion of the extensor carpi radialis tendon
- Circular pad for extra pressure on the muscle in the arm or wrist
- Adjustable compressive strap reduces tension in wrist and finger extensors during gripping Application:
- Tennis elbow
- Golfer's elbow

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Upper limb orthoses

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Elbow Circumference
FBS301335	Oppo - Elbow Brace - Tennis - neoprene - Sml	232 - 257 mm
FBS3013352	Oppo - Elbow Brace - Tennis - neoprene - Med	257 - 283 mm
FBS3013353	Oppo - Elbow Brace - Tennis - neoprene - Lge	283 - 308 mm
FBS3013354	Oppo - Tennis Elbow Brace - neoprene - XL	308 - 333 mm