

## Aspire ComfiMotion Bed Collection - Crook Handle Bed Rail



- Provides multiple grip points to assist repositioning in bed, and safe entry and exit
- Designed to be used in conjunction with ComfiMotion bed collection
- Additional specifications
  - Height of rail from mattress deck: 565 mm
  - BEA830205 width of rail: 680 - 1060 mm
  - BEA830212 width of rail: 685 mm

✓ Contracted DVA product

✓ Contract approved

### Support @ Home

**Function:** Beds and bed equipment

**Category:** Prescribed

Click [here](#) for more information

### Compliance

IEC 60601-2-52:2009

Code	Display Name
BEA830200	Aspire ComfiMotion Care / Relax Bed - Crook Handle Bed Rail
BEA830211	Aspire ComfiMotion Activ Care Bed - Crook Handle Bed Rail