

Aspire ComfiMotion Bed Collection - Crook Handle Bed Rail



- Provides multiple grip points to assist repositioning in bed, and safe entry and exit
- Designed to be used in conjunction with ComfiMotion bed collection
- Additional specifications
 - Height of rail from mattress deck: 565 mm
 - BEA830205 width of rail: 680 - 1060 mm
 - BEA830212 width of rail: 685 mm

✓ Contracted DVA product

✓ Contract approved

Compliance

IEC 60601-2-52:2009



Code	Display Name	RRP	Box Quantity
BEA830200	Aspire ComfiMotion Care / Relax Bed - Crook Handle Bed Rail	167	1
BEA830211	Aspire ComfiMotion Activ Care Bed - Crook Handle Bed Rail	167	1