

Adjustable Night Splint



- Patented Formfit 3D molded technology allows for variable padding thickness in strategic locations for superior comfort and fit
- Soft lycra lined foam interior is plush and breathable, eliminating heat build-up
- Three lycra lined foam padded straps with user friendly centre release buckles ensure immobilisation while maintaining comfort across the flexor tendons
- Adjustable dorsi-flexion straps provide a gentle stretching of the plantar fascia and achilles tendon (adjustable from 90 to 10 degrees dorsi-flexion)

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS301607L	Plantar Fasciitis Splint - Adjustable night splint - Lge