

Rolyan DynaPro Flex Elbow



- Abnormal reflex arc therapy provided by the device assists in neuro-rehabilitation
- Axial hinge with legible goniometer facilitates brace fitting
- Elbow cuff extenders available
- New fabric reduces shear/friction
- Gel elbow pad included

- ✔ Contracted DVA product
- ✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Calf/Arm Circumference	Elbow Circumference
FBS301324	Rolyan - Hinged - DynaPro Flex Elbow - Sml	127 - 203 mm	203 - 254 mm
FBS301324L	Rolyan - Hinged - DynaPro Flex Elbow - Lge	228 - 305 mm	279 - 355 mm
FBS301324M	Rolyan - Hinged - DynaPro Flex Elbow - Med	190 - 254 mm	228 - 292 mm