

K Care Commode Seat - Wedge



- Provides specific pelvic positioning
- Marine Grade vinyl is wrapped, stapled and sealed, making cleaning easier with no seam
- Open and closed configuration available
- Comes with plastic seat runners for easy attachment and removal for cleaning

Warranty Information

- 1 Year



Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_093330084_0103_1_2

Mobile Shower Commode –

Composite: 05_413_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name	Seat Depth	Seat Width
BTC074100	Closed Front Commode Seat - Hand Made - Wedge	430 mm	440 mm