

## **K Care Commode Seat - Wedge**



- Provides specific pelvic positioning
- Marine Grade vinyl is wrapped, stapled and sealed, making cleaning easier with no seam
- Open and closed configuration available
- Comes with plastic seat runners for easy attachment and removal for cleaning

## Warranty Information

• 1 Year

/	$\overline{}$
1	-1
L٦	-
v	٠,

## ndis

**Item Code:** 

05\_093330084\_0103\_1\_2

Mobile Shower Commode –

Composite: 05\_413\_0103\_1\_2

Click **here** for more information

## **Compliance**

AS/NZS 3973:2009

Code	Display Name	Seat Depth	Seat Width
------	--------------	------------	------------

BTC074100 Closed Front Commode Seat - Hand Made - Wedge 430 mm 440 mm