

Darco Softie Post Operative Shoes



- Economical, post-op shoe indicated for soft tissue trauma or surgery.
- Semi-rigid metatarsal shank provides control of forefoot motion to protect after surgery to the metatarsal head area of the forefoot.
- Tri-Laminated upper made of light weight tricot-foam material is breathable and comfortably conforms to bony abnormalities.
- Hook and loop straps are easy to fasten and unfasten with one hand.

- ✓ Contracted DVA product
- ✓ Contract approved

Code	Display Name	Women's Shoe Size	Men's Shoe Size
FBS280230	Darco Softie Post Operative Shoes - Women - Sml	4-6	-
FBS280231	Darco Softie Post Operative Shoes - Women - Med	6.8-8	-
FBS280232	Darco Softie Post Operative Shoes - Women - Lge	8.5-10	-
FBS280233	Darco Softie Post Operative Shoes - Male - Sml	-	6-8
FBS280234	Darco Softie Post Operative Shoes - Male - Med	-	8.5-10
FBS280235	Darco Softie Post Operative Shoes - Male - Lge	-	10.5-12
FBS280236	Darco Softie Post Operative Shoes - Male - XL	-	12.5-14

Code	Display Name	Women's Shoe Size	Men's Shoe Size
------	--------------	-------------------	-----------------