

## Conni Chair Pad - Small



- Fluid Holding Capacity 1000 mls
- Suitable for: vehicles, wheelchairs or general furniture
- Soft polyester top layer, highly absorbent core and breathable PUL waterproof backing
- Fully Washable and re-usable
- Include small rubberised dots in order to prevent the pad from moving



**ndis**

**Item Code:**

03\_092489060\_0103\_1\_1

Click [here](#) for more information

### Warranty Information

- 1 Year

*Conni.*

### Compliance

No applicable compliance standard

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Capacity
CON225020	Chair Pad - Small - Tartan	36	21	480 mm	480 mm	1000 ml
CON225025	Chair Pad - Small - Teal	36	21	480 mm	480 mm	1000 ml
CON255023	Conni Chair Pad - Small - Charcoal -	-	-	480 mm	480 mm	1000 ml