

## Jobst UlcerCare Replacement Liners



- Enables effective management of chronic lower leg wounds
- Provide mild compression on their own (17 mmHg), and 40 mmHg when worn in conjunction with the outer UlcerCare Stocking
- Contains silk for healthy skin, easy donning and patient concordance
- Can be worn 24 hours a day and may be used to manage lower leg oedema during bed rest
- Integrated heel and toe for a very comfortable fit

✓ Contracted DVA product

✓ Contract approved



Code	Display Name	Calf/Arm Circumference	Ankle Circumference
PTS641025	Jobst UlcerCare - Replacement Liners (pkt 3) - 17mmHg - Sml	29-35cm	18-21cm
PTS641030	Jobst UlcerCare - Replacement Liners (pkt 3) - 17mmHg - Med	34-41cm	21-24cm
PTS641040	Jobst UlcerCare - Replacement Liners (pkt 3) - 17mmHg - Lge	40-47cm	24-27cm
PTS641031	Jobst UlcerCare - Replacement Liners (pkt 3) - 17mmHg - XL	46-54cm	27-30cm
PTS641032	Jobst UlcerCare - Replacement Liners (pkt 3) - 17mmHg - XXL	48-56cm	30-33cm