

Askle T1 Universal Microbead Cushion



- The Askle Sante range of positioning cushions allow you to effectively and rapidly secure postures which contribute notably in terms of reducing the pressure placed on body skin and tissue, vulnerable to pressure sores.
- The T1 Universal cushion is suitable for all care situations (to be used alone or with other positioning cushions).
- Suitable for use as a pillow which is compliant with prevention of cross-infection risk.
- Easy to install and positioning by carers.
- Compact for easy storage.
- Ultrasonically welded cover to prevent the penetration of liquids.
- Silver Ag+ IONS – Sanitized Silver fighting the development of micro-organisms.
- PU-Polycarbonate Coated Polyamide fabric for improved performance and durability.
- Oeko-Tex class 2 certified for direct contact with skin.
- Self-extinguishing expanded Polystyrene microbeads for improved safety.
- Vacuum Touch™ concept for a stable distribution of the microbeads within their covers.

