

Transfer Rota Seat



- Reduces the risk of strain or back injury for the carer and user.
- Assists the user in transfers from a seated position.
- Ideal for transfers in a variety of situations such as car, wheelchair and household areas.
- Swivels 360.
- Lightweight and portable.
- Smooth and safe turning action.
- Covered in soft fabric upholstery.
- Nonslip underside.

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Changing body position

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Width	Total Height	Safe Working Load	Unit Weight
TSE672860	Rota Seat	390 mm	50 mm	100 kg	2.4 kg