

## Aspire Bed Moving Trolley



- The Bed Moving Trolley is designed to assist in the safe handling and manoeuvring of one-piece beds (long single to queen sizes) over relatively flat and even terrain.
- It should not be used to carry mattresses or other items at the same time.
- The trolley has a narrow base and easily manoeuvres within small corridors making home bed deliveries safer and easier.
- The trolley requires at least two people to safely handle double or queen size beds.
- More than two people are required when handling any bed over rough or uneven terrain, or when a risk assessment determines that more people are necessary.

 Contract approved

### Support @ Home

**Function:** Beds and bed equipment

**Category:** Prescribed

Click [here](#) for more information



Code	Display Name
BEA824205	Bed Moving Trolley for Single to Queen Size Beds - Adjustable Fastening
BEA824206	Bed Moving Trolley for Single to Queen Size Beds - Winch Upgrade Kit