

Spex Vigour - Lo Active Back Support



New invertible shell design means the Vigour Lo can be used either way up on the wheelchair. Sacral and lumbar height support allowing unrestricted upper body range of movement. Designed for the individual who requires:

- A solid back support with the ability to self-propel
- Improved spinal and pelvic positioning whilst maintaining propelling efficiency
- Minimal postural support
- For clients that have mild torso asymmetries
- Approx. Vigour Lo weight is from 0.8 to 1.5kg

Warranty Information

- 2 Years

✔ Contract approved

Support @ Home

Function: Accessories for wheelchair seating

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Safe Working Load	Backrest Height	Backrest Width
PTB535014	Vigour Kidz - Lo Active Back Support For 10in W/C (6in High) For Low Sacral/Lumbar Support	90 kg	6 in	9 in