

Invacare Birdie 170 Patient Lifter



- Birdie 170 can be folded and unfolded without the need for tools or special skills required. This makes it ideal for the nursing or home care environments where space and time are at a premium
- All attachments on the Birdie 170 are quick release, making it very easy to change when required
- Intelligent control system includes a service reminder, emergency operation, audio battery alarm, and battery indicator
- Available in manual or electric leg
- Safe patient handling with up to 170 kg
- Supplied with yoke attachment
- Electric emergency lowering up and down
- Item Weight: 44 kg

Optional Accessories:

- Pivot Frame (LSS364112)
- Full range of slings

Warranty Information

- 3 Years Electronics
- 2 Years Battery
- 1 Year Handset

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Electric Leg	Safe Working Load	Material
LSS364110	Invacare Birdie Patient Lifter - Electric	1150 mm	680 mm	1740 mm	470 - 1740 mm	580 - 1040 mm	True	180 kg	Steel

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Electric Leg	Safe Working Load	Material
------	--------------	--------------	-------------	--------------	------------	----------------	--------------	-------------------	----------

Leg