

Walking Tutor - Adult



- Provides assistance for those who have difficulty walking.
- Tutors the user into an upright gait.
- Folding design for storage or transportation.
- Height adjustable forearm supports to suit user requirements.
- Adjustable hand grips to suit user forearm length.
- Four swivel castors for manoeuvrability.
- Powder coated mild steel construction.
- Skis or glides available as an option as well as padded arm rest.
- Small size designed for children.

Warranty Information

- 1 Year

 Contract approved

Support @ Home

Function: Walking, manipulated by both arms

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_120606094_0105_1_2

Walking Aids -

Composite:05_416_0105_1_2

Click [here](#) for more information

Compliance

ISO 11199-2:2005

Code	Display Name	Total Length	Total Height	Safe Working Load
WAF708900	Walking Tutor - Adult	515 mm	850 - 1350 mm	100 kg
WAF709000	Walking Tutor - Adult w/ Glides	515 mm	850 - 1350 mm	100 kg
WAF709100	Walking Tutor - Adult w/ Padded Chest Rest	515 mm	850 - 1350 mm	100 kg

Code	Display Name	Total Length	Total Height	Safe Working Load
------	--------------	--------------	--------------	-------------------