

Conni Chair Pad - Large



- Fluid Holding Capacity 1000 mls
- Suitable for: vehicles, wheelchairs or general furniture
- Soft polyester top layer, highly absorbent core and breathable PUL waterproof backing
- Fully Washable and re-usable
- Include small rubberised dots in order to prevent the pad from moving



ndis

Item Code:

03_092489060_0103_1_1

Click [here](#) for more information

Warranty Information

- 1 Year

Conni.

Compliance

No applicable compliance standard

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Capacity
CON225030	Chair Pad - Large - Tartan	43	21	510 mm	610 mm	1000 ml
CON225035	Chair Pad - Large - Teal	43	21	510 mm	610 mm	1000 ml