

TributeWrap Glove



- TributeWrap is an adjustable compression garment with foam technology that nestles into the tissue to promote skin health. It supports easy self-management of lymphoedema and oedema during low activity by comfortably adjusting to your patient's unique shape and lifestyle
- TributeWrap helps maintain gains made during therapy while allowing the flexibility to support further reduction. Simply slide it on, adjust the straps, and cover with Sleep Sleeve for comfortable evening compression

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name
PTA537800	TributeWrap Glove Small Right
PTA537790	TributeWrap Glove Small Left
PTA537780	TributeWrap Glove Medium Right
PTA537770	TributeWrap Glove Medium Left
PTA537760	TributeWrap Glove Large Right
PTA537750	TributeWrap Glove Large Left