

TED - Anti-Embolic - Thigh w/ Belt Stockings



- The TED thigh high with waist belt is very suitable for thin or obese patients, where stocking slippage is often a problem
- One-piece, snap-on belt is easily adjustable
- Pressure break at knee, protecting the critical popliteal area and promoting free leg movement in the recumbent patient
- Available in 5 sizes, with choice of either short, regular or long fit
- Supplied per pair

 Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name	Colour
PTS641425	TED Anti-embolic Stockings - Thigh w/ Belt - Toe Hole - X Sml / Reg	White
PTS641430	TED Anti-embolic Stockings - Thigh w/ Belt - Toe Hole - Sml / Reg	White
PTS641460	TED Anti-embolic Stockings - Thigh w/ Belt - Toe Hole - Med / Reg	White
PTS641470	TED Anti-embolic Stockings - Thigh w/ Belt - Toe Hole - Lge / Reg	White
PTS641500	TED Anti-embolic Stockings - Thigh w/ Belt - Toe Hole - XL / Reg	White
PTS641501	TED Anti-embolic Stockings - Thigh w/ Belt - Toe Hole - X Sml / Long	White
PTS641502	TED Anti-embolic Stockings - Thigh w/ Belt - Toe Hole - Sml / Long	White

Code	Display Name	Colour
PTS641503	TED Anti-embolic Stockings - Thigh w/ Belt - Toe Hole - Med / Long	White
PTS641504	TED Anti-embolic Stockings - Thigh w/ Belt - Toe Hole - Lge / Long	White
PTS641505	TED Anti-embolic Stockings - Thigh w/ Belt - Toe Hole - XL / Long	White