

## Syringe Pump Stand



- 4 point “pig tail” IV bag attachment
- Adjustable pole with pneumatic dampener mechanism to retard adjustment speed
- Stainless steel 2 5mm upright can accept up to 4 syringe pumps / volumatic infusion pumps / PCA pumps
- Push handle
- 4 pin power board AS3760 (earthed)
- Heavy duty orange power cable
- Stable 6 star base (earthed) 600 mm diameter
- 100 mm grey rubber swivel castors

- ✓ Contracted DVA product
- ✓ Contract approved

### Support @ Home

**Function:** Administering medicines

**Category:** Prescribed

Click [here](#) for more information

### Compliance

No applicable compliance standard

## Warranty Information

- 5 Years

Code	Display Name	Total Height
WAT820150	Syringe Pump Stand - Mobile - 6 Star Base + Power Board + Push Handle + 4 Hook	1550 - 2000 mm
WAT820160	Syringe Pump Stand - Mobile - 6 Star Base + 4 Hook - Economy - No P/board or P/handle	1550 - 2000 mm