


LIKO COMFORTSLING PLUS (MOD. 300)



- Comfort Sling Plus™ is a further development of the Comfort Sling™ (Mod. 30) with improved comfort and fit
- Comfort Sling Plus™ enables an upright sitting posture and support for the entire back
- Comfort Sling Plus™ is applied and removed with the patient in a supine position and is then left in the chair after the transfer
- Suits most patients
- Excellent choice when lifting from a flat position
- Gives excellent support for back and head
- Comfortable reclined sitting posture
- Disposable sling intended for single use (do not wash).

Warranty Information

- 1 Year

 Contract approved

ndis

Item Code:

05_123603164_0105_1_2

Hoist –

Composite: 05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394737	ComfortSling Plus 300 - Net Polyester - MEDIUM	200 kg	Net Polyester
LSS394736	ComfortSling Plus 300 - Mesh - XLARGE	500 kg	Mesh
LSS394720	Comfort Sling 350 - Mesh - LARGE	200 kg	Mesh (Loops)
LSS394730	Comfort Sling 350 - Mesh - XLARGE	300 kg	Mesh (Loops)