

Low Back Chair (Height adjustable)



- A versatile, low back orthopaedic chair
- Contoured Backrest – Provides postural support
- Softly padded armrests, seat and back upholstery are comfortable and supportive
- Non-Slip Rubber Tips – Provide stability whilst sitting or transferring
- Easy clean vinyl surfaces
- Height Adjustable – Easy push button adjustment for comfortable positioning and changing needs
- Lightweight Design - lightweight yet sturdy design for easy manoeuvrability



*** Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

| Code | Display Name |
|-------------|------------------------------------|
| H_CHP208000 | Low Back Chair (Height adjustable) |