

## Shoulder Immobiliser Sling



- High quality, durable poly cotton sling with contact closure adjustment on shoulder and body straps.
- Fits right or left arm and helps prevent shoulder rotation.
- Shoulder pad provides additional comfort.
- Ideal for immobilisation and support of the shoulder and elbow joints.

✓ Contract approved

### Support @ Home

**Function:** Upper limb orthoses

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name
FBS280007	Performance Health Shoulder Immobiliser Sling - Medium - Single