

Kerry Swift Lift Homecare Patient Lifter



- Mobile patient lifter for transfer to a seated or lying position
- Premier version of Swift Lift model
- Robust epoxy powder coated frame construction makes it ideal for aged care facilities
- Concave mast provides more leg room for users
- Dual speed hand control operation
- Emergency lowering and raising
- Over-sized push handles for excellent manoeuvrability
- Electric leg gating
- Single hook yoke and sling retaining system gives a smaller lifting footprint
- Suits a large range of slings
- Supplied with battery and charger



ndis

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Electric Leg	Safe Working Load	Material
LSS379750	Kerry Swift Lift Advantage Patient Lifter	1160 mm	670 mm	1240 - 1710 mm	490 - 1710 mm	515 - 870 mm	True	160 kg	Steel