

## JAY - J2 DEEP CONTOUR Backrest



**JAY**<sup>®</sup>

- Designed for the client who has moderate postural needs who needs a lightweight, easy-to-use back support that provides both pelvic and trunk positioning
- Multiple positioning accessories; height adjustability and tilt of up to 15° allow the J2 DC backrest to alter in accordance with the client's changing needs over time
- Deep lateral walls provide excellent lateral trunk stability whilst the shape and size of the backrest minimises upper limb interference and don't compromise seat depth
- All models include: Universal mounting hardware, 2 lateral support wedges, 2 posterior pelvic supports, and Jay AirExchange cover

### Warranty Information

- 2 Years Backrest
- 6 Months Cover

 Contract approved

### Support @ Home

**Function:** Accessories for wheelchair seating

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_181003192\_0105\_1\_2

**Wheelchair Manual – Composite:**

**Wheelchair Powered –**

**Composite:**05\_417\_0105\_1\_2

05\_419\_0105\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	Safe Working Load	Material	Seat Width	Backrest Height	Backrest Width
PTB532220	JAY J2 Deep Contour - 18in	113 kg	Aluminium, Foam	17 - 19 in	19 in	18 in

