

MEDI 3.266 ELASTIC ELBOW SUPPORT



- Highly-elastic, anatomically formed, knitted two-way stretch support
- The compressive effect improves resorption of oedema and effusions
- Rapid reduction of swelling, with pain relief in the affected region

Application:

- Swellings and joint effusions
- Rheumatic and osteoarthritic irritations

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS307080	MEDI 3.266 ELASTIC ELBOW SUPPORT - SIZE I
FBS307090	MEDI 3.266 ELASTIC ELBOW SUPPORT - SIZE II
FBS307100	MEDI 3.266 ELASTIC ELBOW SUPPORT - SIZE III
FBS307110	MEDI 3.266 ELASTIC ELBOW SUPPORT - SIZE IV
FBS307120	MEDI 3.266 ELASTIC ELBOW SUPPORT - SIZE V