

## BC1 High Back Day Chair



- Stylish orthopaedic ergonomic chair
- Padded armrests and padded fully welded one piece shell
- Fully adjustable seat height
- Powder coated steel frame
- Range of upholstered finishes and colours available
- Please specify colour at time of order

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Sitting furniture

**Category:** Under advice

Click [here](#) for more information

### Warranty Information

- 10 Years Structural
- 1 Year Fabric

### NDIS

**Item Code:**

05\_188800051\_0103\_1\_2

Click [here](#) for more information

Code	Display Name	Total Width	Total Height	Safe Working Load	Material	Colour	Seat Depth	Seat Width	Seat To Floor Height
CHP196554	BC1 Highback Day Chair - Fabric - Tweed Charcoal	730 mm	990 - 1180 mm	180 kg	Fabric	Tweed Charcoal	490 mm	510 mm	360 - 510 mm
CHP196556	BC1 Highback Day Chair - Fabric - Tweed Oatmeal	730 mm	990 - 1180 mm	180 kg	Fabric	Tweed Oatmeal	490 mm	510 mm	360 - 510 mm
CHP196558	BC1 Highback Day Chair - Fabric - Tweed	730 mm	990 - 1180	180 kg	Fabric	Tweed Denim	490 mm	510 mm	360 - 510 mm

Code	Display Name	Total Width	Total Height	Safe Working Load	Material	Colour	Seat Depth	Seat Width	Seat To Floor Height
	Denim		mm						
CHP196560	BC1 Highback Day Chair - Fabric - Tweed Stone	730 mm	990 - 1180 mm	180 kg	Fabric	Tweed Stone	490 mm	510 mm	360 - 510 mm
CHP196562	BC1 Highback Day Chair - Vinyl - Ink	730 mm	990 - 1180 mm	180 kg	Vinyl	Vinyl Ink	490 mm	510 mm	360 - 510 mm
CHP196564	BC1 Highback Day Chair - Vinyl - Latte	730 mm	990 - 1180 mm	180 kg	Vinyl	Vinyl Latte	490 mm	510 mm	360 - 510 mm
CHP196568	BC1 Highback Day Chair - Vinyl - Kiwi	730 mm	990 - 1180 mm	180 kg	Vinyl	Vinyl_Kiwi	490 mm	510 mm	360 - 510 mm
CHP196570	BC1 Highback Day Chair - Vinyl - Hollywood	730 mm	990 - 1180 mm	180 kg	Vinyl	Vinyl Hollywood	490 mm	510 mm	360 - 510 mm
CHP196572	BC1 Highback Day Chair - Vinyl - Red	730 mm	990 - 1180 mm	180 kg	Vinyl	Vinyl Red	490 mm	510 mm	360 - 510 mm
CHP196519	BC1 Highback Day Chair - King - Ink Vinyl	-	990 - 1140 mm	220 kg	-	Ink Vinyl	610 mm	610 mm	630 mm
CHP196520	BC1 Highback Day Chair - King - Latte Vinyl	-	990 - 1140 mm	220 kg	-	Latte Vinyl	610 mm	610 mm	630 mm
CHP196525	BC1 Highback Day Chair - Super King - Ink Vinyl	-	990 - 1140 mm	300 kg	-	Ink Vinyl	490 mm	740 mm	630 mm
CHP196526	BC1 Highback Day Chair - Super King - Latte Vinyl	-	990 - 1140 mm	300 kg	-	Latte Vinyl	490 mm	740 mm	630 mm