

Aspire Overbed Table Top - Laminated & Customs



- Stable and secure table top to conveniently accommodate meals, drinks, reading materials, tablet devices and other items within easy arms reach from a bed or chair
- Height adjustment lever is easily accessible to users and carers
- Height fixing knob will lock the height and prevent adjustment by users
- Two locking castors for securing table position by the carer to prevent it rolling away

Warranty Information

- 2 Years



 Contract approved

ndis

Item Code:

05_180315322_0103_1_2

Minor Bed – Composite:

05_428_0103_1_2

Click [here](#) for more information

Compliance

ISO7172

Code	Display Name	RRP	Box Quantity	Total Length	Total Width
BEA006695	Aspire Overbed Table - Laminated & Customs - Top - Sublime Teak	123	1	430 mm	800 mm
BEA006696	Aspire Overbed Table - Laminated & Customs - Top - Milano Walnut	123	1	430 mm	800 mm
BEA006691	Aspire Over Bed Table - Laminated & Customs - Top - Recessed Cream	106	1	430 mm	800 mm
BEA006694	Aspire Over Bed Table - Laminated & Customs - Top - Beech	68	1	430 mm	800 mm
BEA006693	Aspire Overbed Table - Laminated & Customs - Top - AC White	68	1	430 mm	800 mm