

Aspire Signature 2 CUSTOM Lift Recline Chair - Dual Action





Customise your Aspire Signature 2 - Dual Action Lift Recline Chair with upholstery to match your style and needs.

- Healthcare Vinyl The Lustrell Ranger range is high performance vinyl with HealthGuard Antibacterial and Antimicrobial treatment to assist with infection control and hygiene Commercial quality fabric for long lasting durability
- Microsuede Suede like texture combined with high performance characteristics ensures the Eastwood fabric range remains easy to clean and maintain
- Leather La Casa is a luxurious artisan leather with a soft glove-like texture,
 will retain their unique suppleness, provided they are regularly cared for
- **Soft Touch Fabric** The soft touch weave and diverse colour palette makes the Mazza and Manisa ranges a highly practical and functional solution
- · Commercial quality fabric for long lasting durability
- Custom order lead time varies, approximately 5 6 weeks depending on fabric availability
- Custom orders are subject to terms and conditions please see page two on the order form below
- All colour swatches shown may not be a true representation of the actual product and should be used as a guide only

Warranty Information

2 Years

⊘Contract approved



Code	Display Name	Total Width	Total Height	Safe Working Load	Material	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP228280VL	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Large - Healthcare Vinyl (choose colour)	820mm	1110mm	155kg	Vinyl	510- 535mm	540mm	500mm	500mm	650mm
CHP228280MS	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Large - Microsuede (choose colour)	820mm	1110mm	155kg	Microsuede	510- 535mm	540mm	500mm	500mm	650mm
CHP228280LT	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Large - Leather (choose colour)	820mm	1110mm	155kg	Leather	510- 535mm	540mm	500mm	500mm	650mm
CHP228280FA	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Large - Soft Touch Fabric (choose colour)	820mm	1110mm	155kg	Soft Touch Fabric	510- 535mm	540mm	500mm	500mm	650mm
CHP228250VL	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Medium - Healthcare Vinyl (choose colour)	780mm	1050mm	155kg	Vinyl	490- 515mm	510mm	480mm	480mm	620mm
CHP228250MS	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Medium - Microsuede (choose colour)	780mm	1050mm	155kg	Microsuede	490- 515mm	510mm	480mm	480mm	620mm
CHP228250LT	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Medium - Leather (choose colour)	780mm	1050mm	155kg	Leather	490- 515mm	510mm	480mm	480mm	620mm
CHP228250FA	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Medium - Soft Touch Fabric	780mm	1050mm	155kg	Soft Touch Fabric	490- 515mm	510mm	480mm	480mm	620mm



Code	Display Name	Total Width	Total Height	Safe Working Load	Material	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
	(choose colour)									
CHP228220VL	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Small - Healthcare Vinyl (choose colour	740mm	980mm	155kg	Vinyl	460- 385mm	480mm	480mm	480mm	580mm
CHP228220MS	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Small - Microsuede (choose colour)	740mm	980mm	155kg	Microsuede	460- 385mm	480mm	480mm	480mm	580mm
CHP228220LT	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Small - Leather (choose colour)	740mm	980mm	155kg	Leather	460- 385mm	480mm	480mm	480mm	580mm
CHP228220FA	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Small - Soft Touch Fabric (choose colour)	740mm	980mm	155kg	Soft Touch Fabric	460- 385mm	480mm	480mm	480mm	580mm