

## DermaSaver Elbow Tube Protection



- Soft and loose with ample airflow reaching the skin allowing for greater healing ability whilst offering protection
- Cushioned layers of MicroSpring Textile protect fragile skin from friction and pressure that can cause skin tears, bruising and abrasions
- Sizing: simply measure the circumference around the widest part of the forearm. For borderline measurements, choose the larger size to ensure comfort

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Manage tissue integrity

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Calf/Arm Circumference
PTA533304	DermaSaver - Elbow Tube - X Sml	13-18cm
PTA533304S	DermaSaver - Elbow Tube - Sml	18-25cm
PTA533304M	DermaSaver - Elbow Tube - Med	25-30cm
PTA533304L	DermaSaver - Elbow Tube - Lge	30-38cm
PTA533304XL	DermaSaver - Elbow Tube - XL	38-46cm
PTA533305	DermaSaver - Double Elbow Tube - X Sml	13-18cm
PTA533305S	DermaSaver - Double Elbow Tube - Sml	18-25cm

<b>Code</b>	<b>Display Name</b>	<b>Calf/Arm Circumference</b>
PTA533305M	DermaSaver - Double Elbow Tube - Med	25-30cm
PTA533305L	DermaSaver - Double Elbow Tube - Lge	30-38cm
PTA533305XL	DermaSaver - Double Elbow Tube - XL	38-46cm